



NCSSM

trmB Gene Expression in *Haloferax volcanii* in the Presence of Different Carbon Sources

North Carolina School of Science and Mathematics - Courtesy of Schmid Lab at Duke University



Introduction

- A *trmB* deletion mutant strain, $\Delta trmB$, appears to be complemented by glucose in *H. salinarum*. Complementation means the mutation is no longer expressed (Todor, et al.).
- trmB* is a transcription factor protein that activates other genes to assist in gluconeogenesis (the derivation of glucose). This glucose is used to create S-layer cell surface proteins, which *H. salinarum* and *H. volcanii* need to grow
- It is unclear if *H. volcanii* and *H. salinarum* use carbon sources and *trmB* in a similar way (Todor, et al.)

Figure 1.1: *Haloferax volcanii* is a halophilic archaea with a similar genome to previously studied *Halobacterium salinarum*.





Figure 1.2: It is unclear exactly what role *trmB* plays in the gluconeogenesis process of *H. volcanii* (BioRender)

- Determine the effects of different carbon sources on $\Delta trmB$'s ability to generate its own glucose through gluconeogenesis
- It is hypothesized that in the presence of any carbon source, $\Delta trmB$ would exhibit growth similar to wild-type.**
- The optical density of cell cultures is measured to model their growth in the presence of these carbon sources.
- $\Delta trmB$ exhibited the highest level of growth with glucose. Galactose and glycerol were close to the wild-type levels of growth. Xylose had a considerable population increase before tapering off and exhibiting similar patterns as no carbon.

Methods




Figure 2.1: This experiment utilized the inoculation and optical density measurement method used in the Todor, et al. experiment. (BioRender)



Figure 2.2: *Haloferax volcanii* after being spread into single cultures.



Figure 2.3: The bioscreen machine used to collect OD values (Designs)

Results

- In the presence of galactose and glycerol, $\Delta trmB$ grows slightly less than WT, but still with relative success.
- Glucose complemented $\Delta trmB$ the best, leading to growth rates higher than the WT sample
- In the presence of xylose, $\Delta trmB$ grows less than the wild-type, although xylose's initial growth rate was comparable
- Both strains grew poorly in the absence of carbon




Figure 3.1: Growth of $\Delta trmB$ compared to WT, organized by carbon source




Figure 3.2: Ratio between the growth of $\Delta trmB$ and wild-type with indication of standard deviation.

Conclusion

- In the presence of glucose, galactose, and glycerol, the $\Delta trmB$ mutation was complemented.**
- Due to $\Delta trmB$'s gluconeogenesis inhibition, when the samples were grown with no carbon, $\Delta trmB$ grew less than the wild-type.
- With no carbon, the wild-type's only source of glucose is produced through gluconeogenesis. $\Delta trmB$ could not produce any supplemental glucose, which is why it barely grew.




Figure 4.1: Glucose is a more complex compound than xylose, so xylose cannot be derived from glucose. (Jackson)

Figure 4.2: The xylose data suggests *trmB* activates the *XacR* gene and is needed to activate xylose metabolism.

- Further experimentation should be done to discover the differences in catabolizing xylose versus other carbon sources.
 - $\Delta trmB$'s effect on S-layer growth with different carbon sources
 - $\Delta trmB$'s effect on cell morphology
- $\Delta trmB$ affects metabolic processes of one of the oldest living organisms, hinting toward how metabolic pathways evolved.
- Our results support further research into archaeal metabolic pathways.

References

- Andrei AS, et al. Living with salt: metabolic and phylogenetic diversity of archaea inhabiting saline ecosystems. *FEMS Microbiol Lett*. 2012 May;330(1) 1-9. PMID 22339687
- Designs, WebDuck. "Thermo LabSystems Bioscreen C Automated Microbiology Growth Curve Anal." *Used Lab Equipment* | American Laboratory Trading, American Laboratory Training, 2022, <https://americanlaboratorytrading.com/lab-equipment-products/thermo-labsystems-bioscreen-c-automated-microbiology-growth-curve-analysis-system-type-fp-1100-c-15748>.
- Jackson, S., et al. (2002, February 14). Xylose as a nectar carbon source: From biochemistry to ecology. *Comparative Biochemistry and Physiology Part B: Biochemistry and Molecular Biology*. Retrieved January 25, 2022, from <https://www.sciencedirect.com/science/article/pii/S1096495902000283>
- Pohlschroder, Mechtild, et al. *Haloferax volcanii*. *Trends in Microbiology*, Volume 27, Issue 1. 2019. Pages 86-87, ISSN 0966-842X, <https://doi.org/10.1016/j.tim.2018.10.004>.
- Todor, H., et al. (2014). A transcription factor links growth rate and metabolism in the hypersaline adapted archaeon *Halobacterium salinarum*. *Molecular microbiology*, 93(6), 1172-1182.
- Zwietering MH, et al. 1990. Modeling of the bacterial growth curve. *Appl Environ Microbiol*. 1990 Jun;56(6) 1875-81. PMID 16348228

Acknowledgments

- Schmid Lab
- NCSSM
- Duke University

